



**CITY AND COUNTY OF SAN FRANCISCO
PUBLIC HEALTH LABORATORY**
101 Grove Street, Room 419
San Francisco, CA 94102
Tel: (415) 554-2800 Fax: (415) 431-0651
Godfred Masinde, PhD., Lab Director
CLIA ID # 05D0643643

Test Order *GI Panel* PCR

Synonym(s) Biofire GI Panel, Stool NAAT, Stool PCR. This test can detect nucleic acids from the following bacteria, viruses, and parasites from stool specimen from patients with gastrointestinal symptoms.

- Adenovirus F 40/41
- *Campylobacter* spp.
- Cryptosporidium
- *E. coli* O157
- Shiga-like toxin-producing *Escherichia coli* (STEC) stx1/stx2
- Enterotoxigenic *E. coli* (ETEC)
- *Giardia lamblia*
- *Plesiomonas shigelloides*
- *Salmonella* spp.
- *Vibrio* spp.
- *Yersinia enterocolitica*
- Astrovirus
- *Clostridium difficile* toxin A/B
- *Cyclospora cayetanensis*
- *Entamoeba histolytica*
- Enteropathogenic *E. coli* (EPAC)
- Enteroaggregative *E. coli* (EAEC)
- Norovirus GI/GII
- Sapovirus
- Shigella/enteroinvasive *E. coli* (EIEC);
- *Vibrio cholerae*
- Rotavirus A

Methodology	Polymerase chain reaction (PCR)
Acceptable Specimen Type(s) for Testing	Stool specimens should be collected in Cary Blair transport media according to manufacturer's instructions
Transport / Collection Medium	Cary Blair transport media
Storage and Preservation of Specimen	<ul style="list-style-type: none"> • Specimens should be processed and tested with the BioFire GI Panel as soon as possible. Specimen must be refrigerated (2-8 °C) and received by the Laboratory within three days of collection.
Minimum Volume Required	0.2 mL (200 µL)
Additional Collection Instructions	<p>Introduce stool into orange Para-Pak® vial in small amounts such that level does not exceed fill-line on label. Cap securely and shake vial to distribute sample into the Cary-Blair preservative.</p> <p>If possible, allow patient to urinate before collecting stool specimen so as to avoid contaminating the stool specimen with urine. Catch the stool specimen in a clean, empty wide-mouthed container or place plastic wrap over the opening of the toilet bowl to prevent the stool specimen from falling into the bowl.</p> <p>The stool specimen must be placed into the Para-Pak® vial within one hour of the stools production for optimum results. Label the Para-Pak® vial with the patient's name, date of birth, date of collection and time of collection.</p>
Additional Required Information	Requisition is required for all specimen
Send Out?	N/A
Turnaround Time	2 business days from receipt
Testing Restrictions	
Requisition Form(s)	https://www.sfgcdcp.org/wp-content/uploads/2023/04/Lab-Requisition-Form-Updated-3.17.2023-1.pdf https://www.sfgcdcp.org/wp-content/uploads/2023/04/SFPHL-Bacteriology-Requisition-Form-4.10.2023.pdf
Limitations / Notes / Disclaimers	The performance of this test has not been established for patients without signs and symptoms of gastrointestinal illness. Virus, bacteria and parasite nucleic acid may persist in vivo independently of organism viability. Additionally, some organisms may be carried asymptotically. Detection of organism targets does not imply that the corresponding organisms are infectious or are the causative agents for clinical symptoms. The detection of organism nucleic acid is dependent upon proper sample collection, handling, transportation,

storage and preparation. Failure to observe proper procedures in any one of these steps can lead to incorrect results. There is a risk of false-positive and false-negative results caused by improperly collected, transported or and handled specimens.

Updated 5-17-23